



PHILIPS

RESPIRONICS

DreamStation

Auto CPAP

Designed to **empower patients** and facilitate APAP therapy adoption

DreamStation Auto CPAP features

- Clinically validated auto algorithm
- Proactive search algorithm
- Advanced respiratory event detection*
- Auto-Trial
- CPAP-Check
- Opti-Start
- A-Flex
- **NEW** - EZ-Start
- **NEW** - Smart Ramp : 0 to 45 min
- **NEW** - Performance Check
- Mask fit check
- Compliance and waveform reports
- On-board memory
- **NEW** - remote monitoring solution via Bluetooth®
- Oximetry module (optional) to monitor therapy efficacy

Humidification

- Active humidification controlled by the room environment
- Pre-heat humidifier capable
- Heated tube (optional)

DreamStation Auto CPAP benefits

- **Daily Progress Feedback** displays a simple trend of patients' nightly use to help track progress on the path to adherence.
- **EZ-Start** : automatic, personalized adjustments to CPAP pressure can help patients gradually acclimate to their prescribed level of therapy
- **Connectivity** and remote monitoring options including a standard integrated Bluetooth module
- **Performance Check** simplifies in-home device evaluation for easy remote troubleshooting
- **Patient' data management** through detailed compliance and waveform reports

*The Philips Respiration devices detect and log the following events: Open airway apnea (CA), Obstructive airway apnea (OA), Hypopnea (H), respiratory events-related arousal (RERA), flow limitation (FL), vibratory snore (VS), large leak (LL), periodic breathing (PB).

DreamStation Auto CPAP specifications

Part number	Description
INX500S15	DreamStation Auto CPAP with SD card
INX500H15	DreamStation Auto CPAP with humidifier and SD card

Specifications

Therapy modes	Auto, CPAP, CPAP-Check
Therapy features	Auto-Trial, EZ-Start, Opti-Start
Flex pressure relief	0 to 3
Pressure range	4-20 cm H ₂ O
Smart Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Device controls	Color LCD/control wheel/push button
LCD backlight control	Ambient light sensor
Filters	Reusable pollen and disposable ultra-fine
Data storage capacity (minimum)	Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months
Altitude compensation	Automatic
Warranty	2 years
Dimensions	15,7 (L) x 19,3 (W) x 8,4 cm (H) (without humidifier) 29,7 (L) x 19,3 (W) x 8,4 cm (H) (with humidifier)
Weight	1,33 kg (without humidifier), 1,98kg (with humidifier) but including power supply
Power supply	80W/12V/6,67A
Electrical requirements	100-240 VAC, 50/60 Hz, 2.0-1.0 A

Sleep therapy accessories

Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520

Tubing	Part number
Heated tube (15mm)	HT15
Standard tube (15mm)	PR15

Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519

Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module	1120293
DreamStation SpO ₂ module	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499

Connectivity	Part number
DreamStation cellular/SpO ₂ modem	100621C
DreamStation Wifi/SpO ₂ modem	100721W

