

# NIDEK *Medical*

**Proudly Presents**  
The Newest and Most  
Innovative Member  
of the Nuvo Family  
of Oxygen Concentrators



.....Mark 5  
**Nuvo**  
*Lite*

- "Lite" - a mere 30 lbs. (13.6kg)
- Quiet Operation
- Compact
- Lockable Flow Control Valve
- Liter Flow -- 0.125 - 5 LPM

# NIDEK *Medical*



**Metal DISS  
Fitting**



**Lockable Flow  
Control Valve**



**Circuit Breaker  
Reset**

## Specifications

Nuvo Lite Oxygen Concentrator(STD)

**115 Volt - 60 Hz**

520(STD)

**230 Volt - 50 Hz**

920 (STD)

Nuvo Lite Oxygen Concentrator  
With oxygen Monitoring System

525(OCSI)

925 (OCSI)

Average Power Consumption  
Oxygen Concentration ( at 5 LPM)  
Liter Flow  
Sound Level  
Outlet Pressure  
Weight  
Dimensions

300 watts  
90% (+ 5.5 % / -3% )  
0.125 to 5 LPM  
40 dBA typical  
7.0 PSI  
30 lbs  
14"W × 9"D × 23"H

290 watts  
90% (+ 5.5 % / -3% )  
0.125 to 5 LPM  
40 dBA typical  
7.0 PSI  
13.6 Kg  
35 cm × 22 cm × 58 cm

## Accessories and Re-order Items

Compressor Intake Filter  
External Cabinet Filter  
Humidifier  
Cannula with 7' tubing  
Nuvo Lite User's Guide

8400 - 1180  
8400 - 1025  
9012 - 8774  
9012 - 8780  
2010 - 8400

8400 - 1180  
8400 - 1025  
9012 - 8774  
9012 - 8780  
2010 - 8400 CE